


# Centrifuge Tubes



## 8060 Tube, Centrifuge, Conical Bottom, with Rim

Code No.	Specifications		Packing Quantity	
	Capacity	Max.O.D x Length	Shelf pack	Case
8060CTF10	10ml	16.5x105mm	1	100
8060CTF15	15ml	16.5x110mm	1	100




## 8062 Tube, Centrifuge, Conical Bottom, Screw Cap

Code No.	Specifications			Packing Quantity	
	Capacity	Max.O.D x Length	Cap size	Shelf pack	Case
8062CTF10	10ml	16.5x105mm	GCM415-15	10	50
8062CTF30	30ml	28x105mm	GCM410-24	1	10
8062CTF40	40ml	35x105mm	GCM410-28	1	10

## 8064 Tube, Centrifuge, Conical Bottom, Glass Stopper

Code No.	Specifications			Packing Quantity	
	Capacity	Max.O.D x Length	Stopper size	Shelf pack	Case
8064CTF10	10ml	16.5x100mm	14/10	10	50



## 8080 Tube, Centrifuge, Conical Bottom, Graduated


Code No.	Specifications				Packing Quantity	
	Capacity	Max.O.D xLength	Interval	Tolerance	Shelf pack	Case
8080CTF10	10ml	16.5x105mm	0.2ml	±0.2ml	10	50

## 8082 Tube, Centrifuge, Conical Bottom, Screw Cap, Graduated

Code No.	Specifications					Packing Quantity	
	Capacity	Max.O.D xLength	Interval	Tolerance	Cap size	Shelf pack	Case
8082CTF10	10ml	16.5x105mm	0.2ml	±0.2ml	GCM415-15	10	50


## 8084 Tube, Centrifuge, Conical Bottom, Glass Stopper, Graduated

Code No.	Specifications					Packing Quantity	
	Capacity	Max.O.D xLength	Interval	Tolerance	Stopper	Shelf pack	Case
8084CTF10	10ml	16.5x105mm	0.2ml	±0.2ml	14/10	10	50





# Centrifuge Tubes



## 8540 Tube, Centrifuge, Round Bottom, Heavy Duty

Code No.	Specifications		Packing Quantity	
	Capacity	Max.O.D x Length	Shelf pack	Case
8540CTF50	50ml	35x100mm	1	10
8540CTF100	100ml	45x120mm	1	10