



2918 Cover Glass

Code No.	Specifications		Packing Quantity	
	Width x Length	Thickness	Shelf pack	Case
2918COVER18-18	18x18mm	0.13-0.16mm	200	2,000
2918COVER22-22	22x22mm	0.13-0.16mm	100	1,000
2918COVER24-24	24x24mm	0.13-0.16mm	100	1,000
2918COVER18-24	18x24mm	0.13-0.16mm	100	1,000
2918COVER22-24	22x24mm	0.13-0.16mm	100	1,000
2918COVER24-32	24x32mm	0.13-0.16mm	100	1,000
2918COVER24-40	24x40mm	0.13-0.16mm	100	1,000
2918COVER24-45	24x45mm	0.13-0.16mm	100	1,000
2918COVER24-50	24x50mm	0.13-0.16mm	100	1,000
2918COVER24-55	24x55mm	0.13-0.16mm	100	1,000
2918COVER24-60	24x60mm	0.13-0.16mm	100	1,000

- Easy handling ; 50 pieces are packed vertically in a plastic case.
- Shrink packaging prevents dust, hygroscopicity, etc.



2926 Microscope Slide

Code No.	Specifications			Packing Quantity	
	Width x Length	Thickness	Finish	Shelf pack	Case
2926WSLID-P-N	76x26mm	0.9-1.1mm	Polished	100	5,000
2926WSLID-PF-N	76x26mm	0.9-1.1mm	Polished&Frosted	100	5,000

- Pre-cleaned
- 100 pcs are packed in inner shrink pack. Desiccant(Silicage)is included in outer carton.



2937 Microscope Slide

Code No.	Specifications			Packing Quantity	
	Width x Length	Thickness	Finish	Shelf pack	Case
2937WSLID-P	76x26mm	1.1-1.3mm	Polished	100	4,000
2937WSLID-PF	76x26mm	1.1-1.3mm	Polished&Frosted	100	4,000

- Pre-cleaned
- 100 pcs are packed in inner shrink pack. Desiccant(Silicage)is included in outer carton.



2990 Rack for Microscope Slide

Code No.	Specifications		Packing Quantity	
	Width x Length x Height	Slide per Rack	Shelf pack	Case
2990RACK100	172x204x35mm	100 slides		1
2990RACK25	97x83x33mm	25 slides		50



NOVIX- II, Sealing Film

Code No.	Specifications		Packing Quantity	
	Dimensions of Roll		Shelf pack	Case
NOVIX2-10	Width 10cm x Length 40m (Thickness 100 μm)		1	12

- Colorless film improves in sealingness and attachability.
- Best fit to seal containers such as test tubes, beakers, etc.
- Durable at a temperature of -80℃~+66℃.

