

# Disposable Pipette



## 7077 Pipette, Disposable, Serological, Individually Wrapped CORNING

Code No.	CORNING Code No.	Specifications				Packing Quantity	
		Capacity	Graduation Interval	Negative Graduation	Length	Shelf pack (piece/pack)	Case
7077PIPET0.5	7077-5X	0.5ml	0.01ml	0.2ml	215mm	100	400
7077PIPET1-N	7077-1N	1ml	0.01ml	0.2ml	290mm	200	800
7077PIPET2-N	7077-2N	2ml	0.01ml	0.2ml	290mm	180	720
7077PIPET5-N	7077-5N	5ml	0.1ml	1.0ml	290mm	120	720
7077PIPET10-N	7077-10N	10ml	0.1ml	2.0ml	290mm	100	600

- 7077 pipettes are supplied in individually wrapped, sterile and plugged.
- Easy-open wrapping (one side paper, the other side plastic) has color-coded sizes printed directly on the paper for easy identification by size. Colors conform to I.S.O. standards.
- Do not pipet by mouth. We suggest using a mechanical pipetting device.

## 7078 Pipette, Disposable, Serological, Multi-Pack CORNING


Code No.	CORNING Code No.	Specifications				Packing Quantity	
		Capacity	Graduation Interval	Negative Graduation	Length	Shelf pack (piece/bag)	Case
7078PIPET0.5	7078-5X	0.5ml	0.01ml	0.2ml	215mm	25	500
7078PIPET1-0.1-N	7078-1N	1ml	0.1ml	0.2ml	290mm	50	1,000
7078PIPET1-0.01-N	7078-1CN	1ml	0.01ml	0.2ml	290mm	50	1,000
7078PIPET2-N	7078-2N	2ml	0.01ml	0.2ml	290mm	35	700
7078PIPET5-N	7078-5N	5ml	0.1ml	1.0ml	290mm	30	960
7078PIPET10-N	7078-10N	10ml	0.1ml	2.0ml	290mm	20	720

- 7078 pipettes are supplied in multi-packed, sterile and plugged. Designed for large volume user.
- Similar to 7077 pipettes, except multi-packed in plastic bags.
- Do not pipet by mouth. We suggest using a mechanical pipetting device.

## 7079 Pipette, Disposable, Serological, Multi-Pack CORNING

Code No.	CORNING Code No.	Specifications				Packing Quantity	
		Capacity	Graduation Interval	Negative Graduation	Length	Shelf pack (piece/bag)	Case
7079PIPET0.1	7079-1X	0.1ml	0.01ml	0.01ml	230mm	125	500
7079PIPET0.2	7079-2X	0.2ml	0.01ml	0.02ml	230mm	125	500
7079PIPET0.5	7079-5X	0.5ml	0.01ml	0.2ml	215mm	125	500
7079PIPET1-N	7079-1N	1ml	0.01ml	0.2ml	290mm	50	1,000
7079PIPET2-N	7079-2N	2ml	0.01ml	0.2ml	290mm	35	700
7079PIPET5-N	7079-5N	5ml	0.1ml	1.0ml	290mm	30	960
7079PIPET10-N	7079-10N	10ml	0.1ml	2.0ml	290mm	20	720

- 7079 pipettes are supplied in multi-packed, non-plugged and non-sterile.
- Do not pipet by mouth. We suggest using a mechanical pipetting device.






# Disposable Pipette



## 7077B Pipette, Disposable, Shorty, Serological, Individual CORNING

Code No.	CORNING Code No.	Specifications				Packing Quantity	
		Capacity	Graduation Interval	Negative Graduation	Length	Shelf pack (piece /pack)	Case
7077B-PIPET1	7077B-1	1ml	0.01ml	0.5ml	215mm	250	500
7077B-PIPET5	7077B-5	5ml	0.1ml	3.0ml	215mm	200	400
7077B-PIPET10	7077B-10	10ml	0.2ml	4.0ml	215mm	200	400
7077B-PIPET25	7077B-25	25ml	0.2ml	5.0ml	300mm	100	400
7077B-PIPET50	7077B-50	50ml	0.5ml	3.0ml	358mm	25	50

- 7077B pipettes are supplied in individually wrapped, sterile and plugged.
- These shorty pipettes, calibrated "to deliver" with blow-out, are ideal for use in confined areas.
- 5, 10, 25ml sizes have toiled tops and 50ml has a polypropylene top end.
- Do not pipet by mouth. We suggest using a mechanical pipetting device.



## 7078B Pipette, Disposable, Shorty, Serological, Multi-Packed CORNING

Code No.	CORNING Code No.	Specifications				Packing Quantity	
		Capacity	Graduation Interval	Negative Graduation	Length	Shelf pack (piece /pack)	Case
7078B-PIPET1	7078B-1	1ml	0.01ml	0.5ml	215mm	250	500
7078B-PIPET5	7078B-5	5ml	0.1ml	3.0ml	215mm	200	400
7078B-PIPET10	7078B-10	10ml	0.2ml	4.0ml	215mm	200	400
7078B-PIPET25	7078B-25	25ml	0.2ml	5.0ml	300mm	100	400
7078B-PIPET50	7078B-50	50ml	0.5ml	3.0ml	358mm	25	50

- 7078B pipette are supplied in multi-packed, sterile and plugged.
- These shorty pipettes, calibrated "to deliver" with blow-out, are ideal for use in confined areas.
- 5, 10, 25ml sizes have toiled tops and 50ml has a polypropylene top end.
- Do not pipet by mouth. We suggest using a mechanical pipetting device.




## Tube, Disposable, Westergren, Sedimentation Rate, Self-zeroing CHASE Scientific

Code No	Description	Specifications		Packing Quantity	
		I.D.	O.D. x Length	Shelf pack (piece /pack)	Case
CS-830-S0	Tube	2.55mm	5.0x240mm	100	1,000
CS-BLB-16	Aspiration Bulb	-	-	-	1
CS-RACK-15N	Rack	-	-	-	1

- Chase Westergren sedimentation rate tube can be used for either the standard procedure by drawing blood in a liquid sodium citrate vacuum tube or the modified procedure using a saline-diluted EDTA vacuum tube.
- "Self-zeroing, cotton plugged pipette" – No need to set the meniscus – blood cell settle from the bottom of the saturated cotton plug.



# Disposable Pipette



## Pipette, Disposable, Pasteur, Transfer, Plain

Code No.	Specifications				Packing Quantity	
	O.D.	Tip I.D.	Length	Capillary Length	Shelfpack (pieces/inner box)	Case
IK-PAS-5P	7mm	0.9-1.25mm	146mm(5-3/4")	40mm	200	1,000
IK-PAS-9P	7mm	0.9-1.25mm	229mm(9")	115mm	200	1,000

●IWAKI Pasteur type transfer pipettes are excellent for use in blood banks, hematology, chemistry, bacteriology and all general laboratory applications.

●Do not pipet by mouth. We suggest using a mechanical pipetting device.

## Pasteur Pipette Aid

Code No.	Specifications		Packing Quantity	
	Graduation Interval	Maximum volume	ShelfPack (pieces/inner box)	Case
PAS-PIPETER 1.5	0.1ml	1.5ml		1

●This pasteur pipette aid can be sterilized in an autoclave and has chemical resistance to almost all solvents and reagents.