




Vial Rack/Cap



1880 Vial Rack for 1880SV and 1890CV


Code No.	Specification			Packing Quantity	
	Capacity	Hole ID	For Vial	Shelf pack	Case
1880RACK2	60 vials	16mm	1880SV 2/3, 1890 2/3	1	10
1880RACK5	40 vials	20mm	1880SV5/7/10, 1890CV5/7/10	1	10
1880RACK-C	Cover for 1880RACK-2&1880RACK-5			1	10

- Location number is imprinted on partition.
- Dimensions — Body O.D. 189W x 55H x 117D Body I.D. 184W x 51H x 110D
Cover O.D. 193W x 18H x 121D




9996 Cap for 1880SV (Phenol Resin, Rubber Liner)

Code No.	Specification	Packing Quantity	
		Shelf pack	Case
9996CAP13	GCM I 425-13	100	600
9996CAP15	GCM I 425-15	100	600
9996CAP20	GCM I 400-20	50	200
9996CAP24	GCM I 400-24	50	200




9995 Cap for 1890CV (Polypropylene, Polypropylene Form Liner)

Code No.	Specification	Packing Quantity	
		Shelf pack	Case
9995CAP13	GCM I 425-13	100	600
9995CAP15	GCM I 425-15	100	600
9995CAP20	GCM I 400-20	50	200
9995CAP24	GCM I 400-24	50	200
9995CAP33	GCM I 400-33	50	200



9994 Cap for 1870SNV (Polyethylene, Polystyrene Upper Lens)

Code No.	For Vial	Packing Quantity	
		Shelf pack	Case
9994CAP2	1870SNV3,1870SNV5	100	900
9994CAP3	1870SNV10	100	600
9994CAP4	1870SNV15,1870SNV20	50	400
9994CAP5	1870SNV30	50	400
9994CAP7	1870SNV50,1870SNV100	50	200



9997 Cap for 1855SIV (Urea Resin, with Aluminum Coated Cork Packing)

Code No.	Specification	Packing Quantity	
		Shelf pack	Case
9997CAP425-15	GCM I 425-15	100	1,000
9997CAP400-22	GCM I 400-22	100	1,000