




InGenius LHR Gel Doc System IG/LHR-E (high resolution)



Suitable for capturing variety images

Electrophoresis gel (Ethidium Bromide, Coomassie Blue, Silver Stain, SYBR® Gold, SYBR Green, SYBR Safe, GelStar®, Sypro®Red , Sypro Ruby, Sypro Orange, Fluorescein, Rhodamine Red™ , Texas Red®, Pro-Q® Diamond , Deep Purple™)

Agar plates of dark, light or two colour colonies

Cells in flasks

Autoradiographs

DNA, RNA or protein on membranes & etc

Features

Resolution 1.4-5.5 million pixel CCD camera

16 bit images (0-65535 grey levels dynamic range)

Ambient cooling

Filter drawer

UV to visible light converter screen included

Universal darkroom

Fast & automated analysis GeneTools software

20*30 cm transilluminator included

Promotion valid till 30 Jun, 2009

GENE COMPANY LIMITED

TEL: 2896-6283

FAX: 2515-9371

EMAIL: sales@genehk.com