



Epoch *with Take3*

Multi-Volume Spectrophotometer system

Measure samples as low as 2 μ L with the world's most versatile multi-volume, multi-sample Spectrophotometer System



Epoch



Take3 plate

Features:

- Save valuable sample by using **Take3 plate** to quantify up to **16 samples** in **2 μ L** volumes of DNA, RNA and proteins without dilution
- Detection limit down to **2 ng / μ l dsDNA**
- Run 2 μ L samples, or samples contained in a **microplate, standard cuvette**, or the **patented BioCell quartz cuvette**
- Dial in individual wavelengths from **200 - 999 nm**
- Perform rapid qualitative or quantitative analysis in microplates from **6- to 384- wells**
- Measure **kinetic** rates of reaction
- Run **multiplex assays** quickly and easily
- Save space on the benchtop with its very **small footprint**
- Avoid costly maintenance with the robust, thoughtfully engineered optics and hardware
- **Take3 plate** is also compatible to various Biotek microplate reader including **Epoch, PowerwaveXS2** and **Synergy!**

One System, More Choices, More Assays

Gene Company Ltd TEL: 2896-6283 FAX: 2515-9371 EMAIL: sales@genek.com