



**Gene Company Limited**  
 基因有限公司  
 A Gene Group Company



# Diamond Grip Plus™ Powder-Free Latex Examination Gloves


## ● STANDARD EXAMINATION GLOVE



### Features & Benefits:

- Polymer-coated for easy, powder-free donning
- Longer and thicker than typical exam gloves for added protection
- Fully textured for a secure grip in wet or dry conditions
- Premium, low-modulus natural rubber latex for a true fit and reduced hand fatigue
- Reduced protein
- EN 374-2, EN 374-3

EN 374



### Performance and Comfort is a Premium.

**Diamond Grip Plus™** is made from soft, premium natural rubber latex for a comfortable fit. Fully textured for enhanced sensitivity and control. Manufactured with integrated polymer technologies for powder-free donning and product consistency. **Diamond Grip Plus™** exceeds all proposed FDA and EN/ISO standards for exam glove quality and performance.

**GENE COMPANY LIMITED**  
 TEL: 2896-6283  
 FAX: 2515-9371  
 EMAIL: sales@genehk.com

# DIAMOND GRIP PLUS™

POWDER-FREE LATEX EXAMINATION GLOVES

## STANDARD EXAMINATION GLOVE

### GLOVE PROFILE

#### Diamond Grip Plus\*

Length	240mm / 9.4 in.
Cuff Thickness	0.12mm / 4.7 mils
Palm Thickness	0.14mm / 5.5 mils
Finger Thickness	0.18mm / 7.1 mils
Tensile Strength (MPa)	min. 21
Elasticity/Elongation (%)	min. 700
Powder Mass	max. 2 mg / glove
Protein Level	max. 50 µg/g

### GLOVE ATTRIBUTES

Material:	Premium quality natural rubber latex
Type:	Non-sterile
Shape:	Ambidextrous
Color:	White
Sizes:	Extra Small, Small, Medium, Large, Extra Large
Cuff:	Beaded
Glove Interior:	Powder-free
Glove Exterior:	Fully textured
Usage:	For single use only
Packaging:	100 gloves per dispenser box, by weight 10 Boxes per case

\* ASTM D 3578; EN455 I, II, III; EN 374-2, EN 374-3; CE Class I



## Latex Examination Gloves

See for yourself why **Microflex®** gloves are the first choice of exam glove users nationwide. For samples or more information on any of the fine latex and synthetic products from **Microflex®**, contact your **Microflex®** distributor or for the name of a distributor in your area, call:

**775-746-6600** ext. 1726

Or visit our website at:

[www.microflex.com/int](http://www.microflex.com/int)

Fax: 775-746-6586.

PRODUCT CODE	RANGE OF SIZE	DESCRIPTION	CASE QUANTITY
DGP-INT-XS	EXTRA-SMALL	Latex, Powder-free, Examination	1000
DGP-INT-S	SMALL	Latex, Powder-free, Examination	1000
DGP-INT-M	MEDIUM	Latex, Powder-free, Examination	1000
DGP-INT-L	LARGE	Latex, Powder-free, Examination	1000
DGP-INT-XL	EXTRA-LARGE	Latex, Powder-free, Examination	1000

**CAUTION STATEMENT:** This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.



**MICROFLEX®**  
THE MOST TRUSTED NAME IN GLOVES®

[www.microflex.com/int](http://www.microflex.com/int)

©2006 Microflex Corporation  
P.O. Box 32000  
Reno, NV 89533-2000  
All Rights Reserved.  
05.05227B AA 032206