



Gene Company Limited

基因有限公司

A Gene Group Company



SCANSPEED Centrifuges

1524M - High-Speed Microcentrifuge



- Microprocessor controlled operation.
- Automatic alarm system for over-heat and over-speed.
- Automatic door release.
- 5 acceleration and 5 deceleration ramps for sensitive samples.
- Separate "Pulse" key for quick run.
- Program memory for up to 100 programs.

Technical Specifications

	1524M
Max. RPM	15,000 rpm
Max. RCF	21,206 x g
Max. capacity	1.5 ml / 2.0 ml x 24
Run time	≤ 99 min 59 sec or continuous
Acceleration levels	5
Deceleration levels	5
Program memory	100
Power & Frequency	220 V (110 V optional), 50~60 Hz
Dimension (W x D x H) mm	249 x 368 x 240
Weight without rotor	18 Kg

Rotors & Accessories for 1524M



Tube			
Tube capacity (ml)	1.5 / 2.0	0.2	0.5
Adapter	None		
Adapter ϕ x L (mm)	-	6.5 x 23	8 x 31
Max. rpm	15,000		
Max. rcf (g-force)	21,206	16,854	18,363
Radius (mm)	84.3	67	73



Tube		
Tube capacity (ml)	8-tube strip	0.2
Adapter	None	None
Adapter ϕ x L (mm)	-	-
Max. rpm	15,000	
Max. rcf (g-force)	14,187	
Radius (mm)	56.4	

Please contact our sales representative or visit websites, www.genehk.com and www.labogene.com for detailed information.