



Gene Company Limited

基因有限公司

A Gene Group Company



SCANSPEED Centrifuges

Mini - Personal Microcentrifuge



Gene special offer

HK \$ 5,600 / Set

- Compact Design with small footprint.
- Separate "Pulse" key for quick run.
- Rapid acceleration (<12 sec) and deceleration.
- Automatic alarm system for over-heat and over-speed.
- Automatic door release.
- Supplied with an autoclavable fixed-angle rotor and 12 adaptors for 0.2 ml microtubes.
- Optional: a rotor for 4 X 8 tube PCR strips.

Rotors & Accessories for mini



Fixed angle rotor, GAM-1.5-12
12 tubes
Angle : 45°

Tube			
Tube capacity (ml)	1.5 / 2.0	0.2	0.5
Adapter	None		
Adapter ϕ x L (mm)	-	6.5 x 23	8 x 31
Max. rpm	13,500		
Max. rcf (g-force)	12,300	8,863	10,290
Radius (mm)	60.4	43.5	50.5



PCR rotor, GPM-0.2-32
32 tubes(4 rows)
Angle : 45°

Tube		
Tube capacity (ml)	8-tube strip	0.2
Adapter	None	None
Adapter ϕ x L (mm)	-	-
Max. rpm	6,000	
Max. rcf (g-force)	outer row-1,850 / inner row-1,542	
Radius (mm)	outer row-45 / inner row-38.5	

Technical Specifications

	mini
Max. RPM	13,500 rpm
Max. RCF	12,300 x g
Max. capacity	1.5 ml / 2.0 ml x 12
Run time	≤ 30 min
Acceleration time	≤ 12 sec
Deceleration time	≤ 16 sec
Power & Frequency	220 V (110 V optional), 50~60 Hz
Dimension (W x D x H) mm	208 x 245 x 145
Weight without rotor	4.4 Kg

Please contact our sales representative or visit websites, www.genehk.com and www.labogene.com for detailed information

Gene Company Limited Tel: 2896-6283 Fax: 2515-9371 Email: sales@genehk.com